



Dynamic Test Kits for R&D
and Quality Control

FAST FRAP ANTIOXIDANT CAPACITY ASSAY KIT

KF01006

100/200/500 TESTS

96 well plate

BOCKit

A brand of  BioQuoChem

Votre interlocuteur en France, Belgique, Luxembourg et Suisse :
LIBIOS

83, rue Edmond Michelet - 69490 Pontcharra Sur Turdine - France

Tél. : +33 (0)4 74 13 03 02 - Fax : +33 (0)4 74 05 28 25 –

Mail : info@libios.fr - www.libios.fr

INDEX

| | |
|--|----|
| 1. General information | 3 |
| 2. Technical Specifications | 4 |
| 3. Materials..... | 5 |
| 4. Introduction | 7 |
| 5. Assay principle | 8 |
| 6. Sample preparation..... | 9 |
| 7. Assay preparation | 11 |
| 8. Assay protocol..... | 13 |
| 9. Data analysis..... | 14 |
| 10. Interferering substances | 16 |
| 11. Troubleshooting | 17 |
| 12. Researcher notes | 19 |
| 13. Warranties and limitation of liability | 22 |

1. GENERAL INFORMATION

Please read this manual carefully before performing the assay.

PRECAUTIONS

This product is designed for research use only, it is not approved for human or animal use, or clinical diagnosis.

All chemicals should be handled with care and in accordance with laboratory safety practices. Maintain order and cleanliness where dangerous products are used. It is recommended to use basic Personal Protective Equipment. For more information on the risks and preventative measures, check the MSDS available at www.bqckit.com.

Do not use after the expiring date. Store reagents as indicated on the section [Materials](#) on page 5.

TECHNICAL RECOMMENDATIONS

Keep enzymes, heat labile components and samples on ice. Let the components reach room temperature before use.

Invert the bottles a few times to ensure the reagents are well mixed before running the assay. Avoid foaming or bubbles when mixing or reconstituting components. Avoid cross contamination of samples or reagents by changing tips between sample, standard and reagent additions.

Ensure plates are properly sealed or covered during incubation steps. Ensure complete removal of all solutions and buffers from tubes or plates during wash steps. Make sure you have the right type of plate for your detection method of choice. Make sure the heat block/water bath and microplate reader are switched on.

Do not run the standard curve and the samples at different times and do not reuse the calculations of another day. Keep the standard and the samples on the assay for the same amount of time. It is recommended to use a multi-channel pipette if possible.

2. TECHNICAL SPECIFICATIONS



Available sizes:

100 tests: 9 standard, 41 samples

200 tests: 9 standard, 91 samples

500 tests: 9 standard, 241 samples

The calculations are just an estimation assuming that all the samples were tested the same day and that every standard and sample is tested on duplicate. Test number refers to total number of wells to be evaluated.



Volume of sample required:

10 μ l/test



Types of sample compatible:

Biological fluids and beverages.



Linear range:

7.5-250 μ g/ml



Type of detection:

Colorimetric (593 nm)



Sensitivity:

0.01 (OD 593/FRAP μ M)



Time required for the assay:

10 min

3. MATERIALS

MATERIALS SUPPLIED

Keep Reagent A **tightly closed** and **avoid exposure to air**. Never introduce pipette tips in the bottle, instead pour the necessary amount in a beaker container in case not all the tests are performed at once. To ensure stability, BQCKit provides several bottles of Reagent A, so as far as possible, use one bottle at a time. Store kit components as indicated below:

100 tests

| Product | N° bottles | Amount | Storage (before use) | Storage (after use) |
|---------------|------------|--------|----------------------|---------------------|
| Reagent A | 5 bottles | 4.5 ml | RT | RT |
| Standard | 2 vials | Powder | RT | - |
| 96-well plate | - | 1 | - | - |

Each bottle of Reagent A is valid for 20 tests

200 tests

| Product | N° bottles | Amount | Storage (before use) | Storage (after use) |
|---------------|------------|--------|----------------------|---------------------|
| Reagent A | 2 bottles | 22 ml | RT | RT |
| Standard | 2 vials | Powder | RT | - |
| 96-well plate | - | 1 | - | - |

Each bottle of Reagent A is valid for 100 tests

500 tests

| Product | N° bottles | Amount | Storage (before use) | Storage (after use) |
|---------------|------------|--------|----------------------|---------------------|
| Reagent A | 3 bottles | 30 ml | RT | RT |
| Standard | 4 vials | Powder | RT | - |
| 96-well plate | - | 1 | - | - |

Each bottle of Reagent A is valid for 135 tests

MATERIALS NEEDED BUT NOT SUPPLIED

Materials:

- Double distilled water (ddH₂O) as MilliQ
- Pipettes and pipette tips

Instrumentation:

- Microcentrifuge
- Vortex mixer
- Colorimetric microplate reader – equipped with filter for OD 593 nm

4. INTRODUCTION

Antioxidant capacity is an overall ability of organisms or food to catch free radicals and prevent their harmful effect. Antioxidative effect includes protection of cells and cellular structures against the harmful effect of free radicals, especially oxygen and nitrogen. Substances with antioxidative properties are called antioxidants. They are contained in food and food supplements, most commonly in fruits, vegetables, rice, wine, meat, eggs, and another foodstuff of plant and animal origin.


Antioxidative systems include antioxidative enzymes, that is, superoxide dismutase, catalase, glutathione peroxidase, glutathione-S-transferase, and non-enzymatic substrates, such as glutathione, uric acid, lipoic acid, bilirubin, coenzyme Q, vitamin C (L-ascorbic acid), vitamin A (retinol), vitamin E (tocopherol), flavonoids, carotenoids, theine compounds in green tea, and others. Some biomolecules are also considered biologically active and clinically significant antioxidants, for example, transferrin, ferritin, lactoferrin, ceruloplasmin, hemopexin, haptoglobin, and uric acid.

Total antioxidant capacity or TAC has been considered an overall parameter, which alterations have been linked to several conditions as infertility, obesity, cancer and neurodegenerative diseases.

BQC FAST FRAP assay kit is a ready-to-use, easy and highly reproducible assay to test TAC on single antioxidants in aqueous solutions, added to plasma and on beverages with the outstanding feature of being faster than the FRAP classical assay.

5. ASSAY PRINCIPLE

This kit measures the antioxidant activity of compounds that are able to reduce the ferric complex. When the complex is at an acidic pH, in the presence of a suitable antioxidant solution, it is reduced, which shows a maximum of absorbance at 593 nm. This reaction is rapid and proportional to the antioxidant capacity of the sample.



6. SAMPLE PREPARATION

BQCKit have tested the samples indicated below.

| Sample | Preparation required | Dilution factor | Diluent | Long term storage |
|---------------------------|----------------------|-----------------|--------------------|-------------------|
| Plasma | No | - | - | -20 °C |
| RPMI Culture medium + NAC | No | - | - | -20 °C |
| DMEM Culture medium + NAC | No | - | - | -20 °C |
| Gingseng juice | No | - | - | -20 °C |
| Orange juice | No | - | - | -20 °C |
| Pineapple juice | No | - | - | -20 °C |
| Litchi juice | No | 1:2 | ddH ₂ O | -20 °C |
| Apple juice | No | - | - | -20 °C |
| Red berries smoothie | Yes | 1:2 | ddH ₂ O | -20 °C |
| Kiwi smoothie | Yes | - | - | -20 °C |

Samples from abnormal or extreme experimental conditions may require a different dilution factor. For sample preparation instructions refer to the section preparation protocol, on page 9.

Is your sample is not included on this list? Check the [BQCKit Testing Program](#) and get a discount on your next order!

PREPARATION PROTOCOL

Reagents and materials required for sample preparation are not supplied. Take in account the sample volume required per test, refer to section [Technical Specifications](#) on page 4.

Smoothies:



Use a food processor to make the smoothie



Filter through a 0.2 µm membrane filter



Cool on ice to assay or freeze at -20°C

Total time required: 15 min

7. ASSAY PREPARATION

REAGENT PREPARATION

Reagents in this kit are ready to use as supplied.

STANDARD PREPARATION

Add exactly 1 ml of ddH₂O to the standard vials that are going to be used immediately and mix well. Dilute standard 1:10 with ddH₂O. For example: 100 µl standard + 900 µl Reagent ddH₂O.


Prepare the calibration curve in 1.5 mL tubes as shown below.

| | Standard (µl) | ddH ₂ O (µl) | FRAP (µM) |
|---|---------------|-------------------------|-----------|
| 1 | 0 | 100 | 0 |
| 2 | 2.5 | 97.5 | 100 |
| 3 | 5 | 95 | 200 |
| 4 | 7.5 | 92.5 | 300 |
| 5 | 10 | 90 | 400 |
| 6 | 12.5 | 87.5 | 500 |
| 7 | 15 | 85 | 600 |
| 8 | 17.5 | 82.5 | 700 |
| 9 | 20 | 80 | 800 |

Antioxidant activity is expressed as FRAP values (Ferric Reducing Ability of Plasma). These values are related to Fe²⁺ concentration. If preferred, Trolox, ascorbic acid and gallic acid can be used instead, but those are not supplied in this kit.






PLATE SET UP

This scheme is just a recommendation on how to perform the assay. For optimal results, [BQCKit recommends running the standards and the samples at least for duplicate](#), but it is the user's discretion to do so.

|  | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|----|----|----|----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | S1 | S1 | S9 | S9 | C8 | C8 | C16 | C16 | C24 | C24 | C32 | C32 |
| B | S2 | S2 | C1 | C1 | C9 | C9 | C17 | C17 | C25 | C25 | C33 | C33 |
| C | S3 | S3 | C2 | C2 | C10 | C10 | C18 | C18 | C26 | C26 | C34 | C34 |
| D | S4 | S4 | C3 | C3 | C11 | C11 | C19 | C19 | C27 | C27 | C35 | C35 |
| E | S5 | S5 | C4 | C4 | C12 | C12 | C20 | C20 | C28 | C28 | C36 | C36 |
| F | S6 | S6 | C5 | C5 | C13 | C13 | C21 | C21 | C29 | C29 | C37 | C37 |
| G | S7 | S7 | C6 | C6 | C14 | C14 | C22 | C22 | C30 | C30 | C38 | C38 |
| H | S8 | S8 | C7 | C7 | C15 | C15 | C23 | C23 | C31 | C31 | C39 | C39 |

S1-S9: Standard wells, C1-C39: Sample wells


8. ASSAY PROTOCOL

| | | |
|---|---|--|
| 1 |  | Set up the plate design, you can use the BQCKit recommended set up (refer to section Plate set up on page 12) or use your own (refer to section Researcher notes on page 19) |
| 2 |  | Add 10 μ l of the sample or standard previously prepared (refer to sections Sample preparation on page 9 and Standard preparation on page 11). |
| 3 |  | Add 220 μ l of FRAP Working solution previously prepared (refer to section Reagent preparation on page 11) in each sample and standard well. |
| 4 |  | Mix for 4 minutes under continuous stirring |
| 5 |  | Read the absorbance at 593 nm |


9. DATA ANALYSIS

ANALYSIS OF THE STANDARD

If the spectrophotometer or microplate reader was not zeroed with the blank, then average the blank values and subtract the average blank value from the standard and unknown sample values. Create a standard curve by plotting A 593 nm (y-axis) vs. standard, FRAP μM (x-axis).



Below, an example of the calibration curve with the four different antioxidant standards (up to 125 μM) that can be used with this kit can be seen.



ANALYSIS OF THE SAMPLE

Determine the unknown sample concentration using the standard curve from the assayed sample value. Average the OD for the replicates and then apply:


$$\text{FRAP } \mu\text{M} = \left(\frac{\text{OD 593 nm-intercept}}{\text{slope}} \right) * \text{dilution factor}$$

Usual values obtained on the samples:

| Sample type | Range of values (FRAP μM) |
|----------------------------|---------------------------------------|
| Plasma | 500-2000 |
| Fruit Juices and Smoothies | 400-1500 |

Those values are merely informative and can be affected by multiple factors.h

Some results obtained by BQCKit are shown below:



10. INTERFERING SUBSTANCES

The following substances have been found to interfere with the assay:


- Non-antioxidant reducing substances
- Strong basic substances.


11. TROUBLESHOOTING


| Problem | Cause | Solution |
|--|--|---|
| Assay not working | Use of ice-cold buffer | Buffers must be at room temperature |
| | Plate read at incorrect wavelength | Check the wavelength and filter settings of the instrument |
| | Use of a different 96 well-plate | Colorimetric: Clear plates, Fluorometric: black wells/clear bottom plate |
| Sample with erratic readings | Samples not deproteinized (if indicated on protocol) | Use TCA precipitation protocol for deproteinization |
| | Cells/Tissue samples not homogenized completely | Use Dounce homogenizer, increase number of strokes |
| | Samples used after multiple free/thaw cycles | Aliquot and freeze samples if needed to use multiple times |
| | Use of old or inappropriately stored samples | Use fresh samples or store at -80°C (after snap freeze in liquid nitrogen) till use |
| Lower/Higher readings in samples and standards | Presence of interfering substances in the sample | Check protocol for interfering substances |
| | Improperly thawed components | Thaw all components completely and mix gently before use |
| | Allowing reagents to sit for extended times on ice | Always thaw and prepare fresh reaction mix before use |

| | | |
|--|--|--|
| | Incorrect incubation times or temperatures | Verify correct incubation times and temperatures in protocol |
| Standard readings do not follow a linear pattern | Pipetting errors in standard or reaction mix | Avoid pipetting small volumes (<5 µl) and prepare a master mix whenever possible |
| | Air bubbles formed in well | Pipette gently against the wall of the tubes |
| | Standard stock is at incorrect concentration | Always refer to dilutions on the protocol |
| Unanticipated results | Measured at incorrect wavelength | Check equipment and filter setting |
| | Samples contain interfering substances | Troubleshoot if it interferes with the kit |
| | Sample readings above/below the linear range | Concentrate/Dilute sample so it is within the linear range |

12. RESEARCHER NOTES

|  | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A | | | | | | | | | | | | |
| B | | | | | | | | | | | | |
| C | | | | | | | | | | | | |
| D | | | | | | | | | | | | |
| E | | | | | | | | | | | | |
| F | | | | | | | | | | | | |
| G | | | | | | | | | | | | |
| H | | | | | | | | | | | | |

|  | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A | | | | | | | | | | | | |
| B | | | | | | | | | | | | |
| C | | | | | | | | | | | | |
| D | | | | | | | | | | | | |
| E | | | | | | | | | | | | |
| F | | | | | | | | | | | | |
| G | | | | | | | | | | | | |
| H | | | | | | | | | | | | |

|  | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|
| A | | | | | | | | | | | | |
| B | | | | | | | | | | | | |
| C | | | | | | | | | | | | |
| D | | | | | | | | | | | | |
| E | | | | | | | | | | | | |
| F | | | | | | | | | | | | |
| G | | | | | | | | | | | | |
| H | | | | | | | | | | | | |

13. WARRANTIES AND LIMITATION OF LIABILITY

Our partner Bioquochem shall not in any event be liable for incidental, consequential or special damages of any kind resulting from any use or failure of the products, even if Bioquochem has been advised of the possibility of such damage including, without limitation, liability for loss of use, loss of work in progress, downtime, loss of revenue or profits, failure to realize savings, loss of products of buyer or other use or any liability of buyer to a third party on account of such loss, or for any labor or any other expense, damage or loss occasioned by such product including personal injury or property damage is caused by Bioquochem's gross negligence. Any and all liability of Bioquochem hereunder shall be limited to the amounts paid by the buyer for the product.

Buyer's exclusive remedy and Bioquochem's sole liability hereunder shall be limited to a refund of the purchase price, or the replacement of all material that does not meet our specifications.

Said refund or replacement is conditioned on buyer giving written notice to Bioquochem within 30 days after the arrival of the material at its destination.

Expiration date: 3 months from the date of delivery

Votre interlocuteur en France, Belgique, Luxembourg et Suisse :
LIBIOS

83, rue Edmond Michelet - 69490 Pontcharra Sur Turdine - France

Tél. : +33 (0)4 74 13 03 02 - Fax : +33 (0)4 74 05 28 25 –

Mail : info@libios.fr - www.libios.fr